



Timestrip Blood Temp 6 Indicators

- Easy to read and use
- Inert until activation
- Activation/Breach Indication
- Serialisation

How long did a temperature breach last during transport or storage of the blood bag? Knowing how long blood products have been kept at elevated temperatures is vital to hospital and blood bank procedures.

Timestrip Blood Temp 6 irreversible storage temperature indicators will tell you how long a blood product has been exposed to higher temperatures.



Ease of use.

Timestrip Blood Temp 6 is inert prior to activation. They can be shipped and stored at ambient temperatures and have no special storage requirements.

Blood Temp 6 has two windows: an activation window that is yellow prior to activation and a breach window that is white prior to activation. Field activate at room temperature. Activation is possible only when the surrounding temperature is above the stated threshold temperature.

Clear activation.

Firmly squeeze the blister located on top of the indicator. Once activated, a green activation line appears in the yellow activation window and continues to fill the window turns completely green. The breach window will remain white. The Timestrip Blood Temp 6 can be directly applied by removing the adhesive liner on the back. Accurate to +/- 0.5°C.

Communication.

Once activated, Blood Temp 6 will continuously monitor temperature. When the core temperature of the blood bag breaches 6°C for Blood Temp 6 an immediate indication of temperature abuse is recorded: indication of arming window will remain green, while the breach window will turn blue.

The breach window is irreversible and cannot be reset or changed.



Inactive



Active – threshold
not exceeded



Active – threshold
exceeded