





# Forehead Thermometers



- **FDA registered**
- **Liquid Crystal temperature monitoring**
- **Adhesive-backed paper or mylar forms**
- **Single or multi-use**
- **Ideal for children, medical centres, or quarantined patients**

During any illness outbreak, early fever detection not only helps prevent the spread of illnesses, it also ensures faster treatment. Liquid crystal technology-based, FDA registered forehead thermometers provide a simple, cost-effective tool for monitoring body temperature.



## **Prepare for the unexpected**

Forehead thermometers offer a fast, accurate way to monitor core body temperatures and take prompt preventive action. You can obtain an accurate temperature in approximately 15 seconds. It's simply a matter of watching the liquid crystal display until the colors stop changing. And since you can take temperatures non-invasively,

the forehead thermometers are ideal for use with children, medical centres or quarantined patients.

## **Accurate and effective temperature readings**

Backed with years of proven research, manufacturing support and designed for extended shelf life, these forehead thermometers are essential elements in any pandemic, flu or emergency first aid kit.



## Specifications

	C-8702	C-8786
<b>Application</b>	Reusable, non-adhesive back	Single-use, adhesive back
<b>Format</b>	Dual-scale	Dual-scale
<b>Operating Temperature Range</b>	35°C to 40°C / 95°F to 104°F	35°C to 40°C / 95°F to 104°F
<b>Dimensions (LxW)</b>	74.9mm x 11.2mm	44.45mm x 12.7mm
<b>Accuracy</b>	±1°C / 1.8°F	±1°C / 1.8°F
<b>Classification</b>	Class II Medical Device	Class II Medical Device
<b>Shelf Life</b>	Extended	Extended
<b>Configuration</b>	Standard or Customisable	Standard or Customisable

## Directions for Use

Thermometers should be used indoor and at room temperature (approximately 20°C - 30°C / 68° - 86°F).

Avoid direct sunlight or strong heat when taking temperature. Do not take temperature if individual has been eating, drinking or exercising within the last 30 minutes due to the possibility of temperature fluctuation.

1. If using the disposable thermometer, remove backing.
2. Hold thermometer at both ends, pressing firmly on dry forehead for approximately 15 seconds.
3. Read thermometer while on forehead  
**Green:** correct temperature

Note: If you do not see green, blue or tan colours will appear:

**Blue:** add 0.5°C / **Tan:** subtract 0.5°C

4. The reusable thermometer may be cleaned with a soft cloth.

